



Create Mac Icon to Start TG Program


1. Select Search icon, and type "automator"




2. Select "Application" and click "Choose"



3. Type "run" to shorten list and drag "Run Shell Script" onto work area




4. Select “as arguments”, then enter text as show, and finally select File->Save



The text lines entered are below. Note there is a space between “/anaconda3/bin/pythonw” and “/anaconda3/Photometry/TG_V6.4.pyw”

```
cd /anaconda3/Photometry  
/anaconda3/bin/pythonw /anaconda3/Photometry/TG_V6.4.pyw  
exit 0
```

5. Type name for icon, be certain “Applications” is selected, the click Save



6. Go to Application folder (or click on LaunchPad) – and find the new icon. Click the new icon to start the program. NOTE: THE PROGRAM MAY BE SLOW TO OPEN – 15-30 SECONDS OR MORE. BE PATIENT.

