



Create Mac Icon to Start TG Program


1. Select Search icon, and type "automator"



2. Select "Application" and click "Choose"



3. Type "run" to shorten list and drag "Run Shell Script" onto work area



4. Select “as arguments”, then enter text as shown, and finally select File->Save

Text to enter:

```
cd /users/yourusername/opt/anaconda3/photometry  
/Users/yourusername/opt/anaconda3/bin/pythonw /users/yourusername/opt/anaconda3/photometry/TG_V6.8.pyw  
Exit 0
```

Note: Leave a space between ‘pythonw’ and ‘/users’

